

# NextSling™

## THE NEXT-GENERATION SHOULDER SLING

### Left Shoulder Stickers Instructions


The NextSling™ comes fitted for the Right Shoulder so placing stickers on the back of the shoulder buckles allow it to be easily switched for the left shoulder. Place an Orange Sticker on the back of the shoulder buckle next to the right elbow and a Blue sticker on the back of the shoulder buckle next to the right wrist. The stickers are on the other side.





### Switch from Right Shoulder to Left Shoulder


Unbuckle both shoulder straps and both waist straps. Turn the arm brace around so that it fits the left shoulder. The sticker color on the shoulder buckle will match the color of the buckle being inserted into it. Buckle the shoulder strap with the Orange buckle to the shoulder buckle with the Orange sticker. Buckle the shoulder strap with the Blue buckle to the shoulder buckle with the Blue sticker. Buckle both waist straps into place. The Orange Buckle should now be next to the Left Shoulder Biceps Cuff.


### Put on Left Shoulder Instructions


- 


1 Unbuckle left arm shoulder strap and left arm waist strap.
- 


2 Place right shoulder strap over right shoulder and loosen shoulder strap velcro to desired length (adjust as necessary).
- 


3 Buckle left waist strap.
- 


4 Remove velcro on left biceps cuff.
- 


5 Remove left buckle velcro on left forearm.
- 

6 Put left arm into sling (adjust right shoulder strap as necessary for ease of application).
- 

7 Velcro left buckle to left forearm.
- 

8 Velcro left biceps cuff.
- 

9 Put left shoulder strap over left shoulder and buckle.
- 

10 Adjust shoulder strap length as necessary for snug fit. The shoulder strap can be cut to the desired length and the velcro piece can be removed and reattached.
- 

11 Adjust back pad fit as necessary for snug fit.

Remove Shoulder Sling by Reversing the Above Steps.