

# NextSling™

## THE NEXT-GENERATION SHOULDER SLING

### Right Shoulder Stickers Instructions

Place a Blue Sticker on the front of the shoulder buckle next to the right elbow and an Orange sticker on the front of the shoulder buckle next to the right wrist. Only place stickers on the front of the shoulder buckles. The sticker color on the shoulder buckle will match the color of the buckle being inserted into it to make it easy to tell how the sling wears.



### Put on Right Shoulder Sling Instructions

- 

1 Unbuckle right arm shoulder strap and right arm waist strap.


- 

2 Place left shoulder strap over left shoulder and loosen left shoulder strap velcro to desired length (adjust as necessary).


- 

3 Buckle right waist strap.


- 

4 Remove velcro on right biceps cuff.


- 

5 Remove right buckle velcro on right forearm.


- 

6 Put right arm into sling (adjust left shoulder strap as necessary for ease of application).


- 

7 Velcro right buckle to right forearm.


- 

8 Velcro right biceps cuff.


- 

9 Put right shoulder strap over right shoulder and buckle.


- 

10 Adjust shoulder strap length as necessary for snug fit. The shoulder strap can be cut to the desired length and the velcro piece can be removed and reattached.
- 

11 Adjust back pad fit as necessary for snug fit.

Remove Shoulder Sling by Reversing the Above Steps.

